

Prasugrel metabolite 1

Chemical Properties

CAS No. :	239466-86-5
Formula:	C ₂₁ H ₂₅ FN ₂ O ₅ S ₂
Molecular Weight:	468.56
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Prasugrel metabolite 1 (M7 (R-119251)) is one of the inactive metabolites of Prasugrel in the human body.
Targets(IC50)	Drug Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1342 mL	10.671 mL	21.342 mL
5 mM	0.4268 mL	2.1342 mL	4.2684 mL
10 mM	0.2134 mL	1.0671 mL	2.1342 mL
50 mM	0.0427 mL	0.2134 mL	0.4268 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481