


## tert-Butyl N-methyl-N-(piperidin-4-ylmethyl)carbamate

## Chemical Properties

CAS No. :	138022-04-5
Formula:	C <sub>12</sub> H <sub>24</sub> N <sub>2</sub> O <sub>2</sub>
Molecular Weight:	228.34
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	tert-Butyl N-methyl-N-(piperidin-4-ylmethyl)carbamate is a PROTAC linker employed in the synthesis of PROTAC molecules.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3794 mL	21.8972 mL	43.7943 mL
5 mM	0.8759 mL	4.3794 mL	8.7589 mL
10 mM	0.4379 mL	2.1897 mL	4.3794 mL
50 mM	0.0876 mL	0.4379 mL	0.8759 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481