

Fmoc-MeVal-OSu

Chemical Properties

CAS No. :	863971-49-7
Formula:	C ₂₅ H ₂₆ N ₂ O ₆
Molecular Weight:	450.49
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Fmoc-Val-Ala-Gly is an ADC Linker utilized in the synthesis of antibody-drug conjugates (ADC).
Targets(IC50)	ADC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2198 mL	11.099 mL	22.1981 mL
5 mM	0.444 mL	2.2198 mL	4.4396 mL
10 mM	0.222 mL	1.1099 mL	2.2198 mL
50 mM	0.0444 mL	0.222 mL	0.444 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481