

NSD2-IN-5

Chemical Properties

CAS No. : 2604513-53-1
Formula: C₂₇H₂₅F₂N₉
Molecular Weight: 513.56
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	NSD2-IN-5 (Compound 198) is an NSD2 inhibitor with an IC ₅₀ value ranging from 0.001 μM to 0.01 μM. It is utilized for cancer research.
Targets(IC ₅₀)	Histone Methyltransferase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9472 mL	9.736 mL	19.4719 mL
5 mM	0.3894 mL	1.9472 mL	3.8944 mL
10 mM	0.1947 mL	0.9736 mL	1.9472 mL
50 mM	0.0389 mL	0.1947 mL	0.3894 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481