

Methoctramine

Chemical Properties

CAS No. :	104807-40-1
Formula:	C36H62N4O2
Molecular Weight:	582.92
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Methoctramine is a potent, heart-selective M2 muscarinic receptor antagonist. In vivo, Methoctramine can inhibit muscarine-induced bradycardia. This compound is applicable in cardiovascular disease research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7155 mL	8.5775 mL	17.155 mL
5 mM	0.3431 mL	1.7155 mL	3.431 mL
10 mM	0.1716 mL	0.8578 mL	1.7155 mL
50 mM	0.0343 mL	0.1716 mL	0.3431 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481