

CY7.5-Alginate

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	CY7.5-Alginate is a CY7.5-labeled sodium alginate. This compound merges the polysaccharide characteristics of sodium alginate with the fluorescent properties of CY7.5, playing a crucial role in fields such as fluorescence labeling, cell imaging, and drug delivery research (Ex/Em = 770/820 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481