


4-[[[(1,1-Dimethylethoxy)carbonyl]amino]benzeneacetic acid

Chemical Properties

CAS No. :	81196-09-0
Formula:	C ₁₃ H ₁₇ NO ₄
Molecular Weight:	251.28
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	4-[[[(1,1-Dimethylethoxy)carbonyl]amino]benzeneacetic acid is a type of PROTAC linker employed in the synthesis of PROTAC molecules.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9796 mL	19.8981 mL	39.7962 mL
5 mM	0.7959 mL	3.9796 mL	7.9592 mL
10 mM	0.398 mL	1.9898 mL	3.9796 mL
50 mM	0.0796 mL	0.398 mL	0.7959 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481