

GRL-03022

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	GRL-03022 is a potent and highly selective non-covalent inhibitor of G protein-coupled receptor kinase 5 (GRK5), with an IC ₅₀ of 60 nM. It demonstrates over 16,000-fold selectivity for GRK5 compared to GRK2 (IC ₅₀ > 1000 μM). GRL-03022 is applicable in research related to heart failure and cancer.
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481