

## STING agonist-49-PAB-Ala-Val-CO-C2-mal

## Chemical Properties

CAS No. :	2803608-03-7
Formula:	C62H73N17O13
Molecular Weight:	1264.37
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	STING agonist-49-PAB-Ala-Val-CO-C2-mal is a peptide-cSTING linker-payload, used as a Drug-Linker Conjugate for ADC (Antibody-Drug Conjugate). It serves as a component for ADC synthesis.
Targets(IC50)	Drug-Linker Conjugates for ADC

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7909 mL	3.9545 mL	7.9091 mL
5 mM	0.1582 mL	0.7909 mL	1.5818 mL
10 mM	0.0791 mL	0.3955 mL	0.7909 mL
50 mM	0.0158 mL	0.0791 mL	0.1582 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481