

SOD1-aggregation-IN-1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SOD1-aggregation-IN-1 (compound 20) is a potent inhibitor of SOD aggregation. It effectively inhibits the aggregation of mutant SOD1 induced by MG132 in PC12-SOD1G85RYFP cells, with an EC50 value of 2.63 μ M. SOD1-aggregation-IN-1 is applicable for research in amyotrophic lateral sclerosis (ALS).
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481