

B-Raf-IN-20

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	B-RafIN 20 (Compound Z-001) is a novel B-RAF inhibitor with an IC50 value of 37.80 nM. It binds to the DFG-out/C-helix-in conformation of the B-RAF kinase and exhibits anticancer activity against non-small cell lung cancer.
Targets(IC50)	Raf

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481