

Asparaginate calcium

Chemical Properties

CAS No. :	39162-75-9
Formula:	C ₄ H ₆ NO ₄ .1/2Ca
Molecular Weight:	152.14
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Asparaginate calcium is a calcium supplement that can enhance serum calcium levels.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.5729 mL	32.8645 mL	65.7289 mL
5 mM	1.3146 mL	6.5729 mL	13.1458 mL
10 mM	0.6573 mL	3.2864 mL	6.5729 mL
50 mM	0.1315 mL	0.6573 mL	1.3146 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481