

Human Serum Albumin-ICG

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Human Serum Albumin-ICG (HSA-ICG) is a serum albumin labeled with Indocyanine green (ICG). It plays a crucial role in areas such as fluorescence labeling and drug delivery research (Ex/Em = 785/813 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481