

## Tyrosine ethyl ester

## Chemical Properties

CAS No. :	34081-17-9
Formula:	C <sub>11</sub> H <sub>15</sub> NO <sub>3</sub>
Molecular Weight:	209.25
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Tyrosine ethyl ester exhibits antioxidant activity.
Targets(IC50)	ROS Kinase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.779 mL	23.8949 mL	47.7897 mL
5 mM	0.9558 mL	4.779 mL	9.5579 mL
10 mM	0.4779 mL	2.3895 mL	4.779 mL
50 mM	0.0956 mL	0.4779 mL	0.9558 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481