

TGR5 agonist 8

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

TGR5 agonist 8 (Compound 7i) is a TGR5 activator with an EC₅₀ of 571 nM. It significantly reduces blood glucose levels and food intake in mice with acute diabetes and suppresses appetite over a longer period. TGR5 agonist 8 is applicable for research in diabetes and obesity.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481