


tert-Butyl 4-(4-(hydroxymethyl)phenoxy)piperidine-1-carboxylate

Chemical Properties

CAS No. :	321337-38-6
Formula:	C17H25NO4
Molecular Weight:	307.39
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	tert-Butyl 4-(4-(hydroxymethyl)phenoxy)piperidine-1-carboxylate is a PROTAC linker employed in the synthesis of PROTAC molecules.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2532 mL	16.266 mL	32.532 mL
5 mM	0.6506 mL	3.2532 mL	6.5064 mL
10 mM	0.3253 mL	1.6266 mL	3.2532 mL
50 mM	0.0651 mL	0.3253 mL	0.6506 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481