

## KQSSKALQR

## Chemical Properties

CAS No. :	164031-33-8
Formula:	C43H80N16O14
Molecular Weight:	1045.21
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	KQSSKALQR is a BCR-ABL peptide that can induce a CD8+ T cell response and is applicable in the research of chronic myeloid leukemia (CML) vaccines.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9567 mL	4.7837 mL	9.5675 mL
5 mM	0.1913 mL	0.9567 mL	1.9135 mL
10 mM	0.0957 mL	0.4784 mL	0.9567 mL
50 mM	0.0191 mL	0.0957 mL	0.1913 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481