

Iodipamide meglumine

Chemical Properties

CAS No. :	3521-84-4
Formula:	C ₃₄ H ₄₈ I ₆ N ₄ O ₁₆
Molecular Weight:	1530.2
Storage:	Keep away from direct sunlight Store at -20°C <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Iodipamide (Adiopdone) meglumine is an iodine-containing contrast agent primarily used for X-ray imaging, specifically in cholecystography and cholangiography.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6535 mL	3.2675 mL	6.5351 mL
5 mM	0.1307 mL	0.6535 mL	1.307 mL
10 mM	0.0654 mL	0.3268 mL	0.6535 mL
50 mM	0.0131 mL	0.0654 mL	0.1307 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481