


rel-(1S,2R)-2-(((tert-butoxycarbonyl)amino)methyl)cyclopropane-1-carboxylic acid

Chemical Properties

CAS No. :	1638768-80-5
Formula:	C10H17NO4
Molecular Weight:	215.25
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	rel-(1S,2R)-2-(((tert-butoxycarbonyl)amino)methyl)cyclopropane-1-carboxylic acid is a PROTAC linker utilized in the synthesis of PROTAC molecules.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6458 mL	23.2288 mL	46.4576 mL
5 mM	0.9292 mL	4.6458 mL	9.2915 mL
10 mM	0.4646 mL	2.3229 mL	4.6458 mL
50 mM	0.0929 mL	0.4646 mL	0.9292 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481