

Exendin-4, Cy5 labeled

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Exendin-4, Cy5 labeled (Exendin-4-Cys(Cy5)) is an Exendin-4 linked covalently with a Cy5 fluorophore, functioning as a GLP-1 receptor agonist. This compound enables in vivo imaging of β cells, offering a means to assess the dynamic expression of GLP-1R, particularly in type 2 diabetes models.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481