


## 1-Boc-4-(4-Carboxybenzyl)piperazine

## Chemical Properties

CAS No. :	479353-63-4
Formula:	C <sub>17</sub> H <sub>24</sub> N <sub>2</sub> O <sub>4</sub>
Molecular Weight:	320.39
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	1-Boc-4-(4-Carboxybenzyl)piperazine is a PROTAC linker utilized in the synthesis of PROTAC molecules.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1212 mL	15.606 mL	31.212 mL
5 mM	0.6242 mL	3.1212 mL	6.2424 mL
10 mM	0.3121 mL	1.5606 mL	3.1212 mL
50 mM	0.0624 mL	0.3121 mL	0.6242 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481