

HDAC6 ligand-Linker Conjugate 3

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	HDAC6ligand-Linker Conjugate 3 consists of an HDAC6 ligand and a PROTAC linker. It is used in the synthesis of the PROTACHDAC6 degrader 7.
Targets(IC50)	Target Protein Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481