

Methyl 15-bromopentadecanoate

Chemical Properties

CAS No. :	41240-56-6
Formula:	C16H31BrO2
Molecular Weight:	335.32
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Methyl 15-bromopentadecanoate is a PROTAC linker utilized in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9822 mL	14.9111 mL	29.8223 mL
5 mM	0.5964 mL	2.9822 mL	5.9645 mL
10 mM	0.2982 mL	1.4911 mL	2.9822 mL
50 mM	0.0596 mL	0.2982 mL	0.5964 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481