

## CRBN ligand-701

## Chemical Properties

CAS No. :	2925072-72-4
Formula:	C14H16N4O3
Molecular Weight:	288.31
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	CRBN ligand-701 is an E3 ubiquitin ligase cereblon (CRBN) ligand designed for recruiting cereblon proteins. It can be linked via a linker to a target protein ligand to form a PROTAC.
Targets(IC50)	Ligands for E3 Ligase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4685 mL	17.3424 mL	34.6849 mL
5 mM	0.6937 mL	3.4685 mL	6.937 mL
10 mM	0.3468 mL	1.7342 mL	3.4685 mL
50 mM	0.0694 mL	0.3468 mL	0.6937 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481