

## 6-Iodohept-1-yne

## Chemical Properties

CAS No. :	2468-56-6
Formula:	C <sub>6</sub> H <sub>9</sub> I
Molecular Weight:	208.04
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	6-Iodohept-1-yne is a PROTAC linker employed in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8068 mL	24.0338 mL	48.0677 mL
5 mM	0.9614 mL	4.8068 mL	9.6135 mL
10 mM	0.4807 mL	2.4034 mL	4.8068 mL
50 mM	0.0961 mL	0.4807 mL	0.9614 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481