

T-02

## Chemical Properties

CAS No. : 26115-72-0  
Formula: C7H15NO3Si  
Molecular Weight: 189.29  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	T-02 is an isocyanate silane coupling agent.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2829 mL	26.4145 mL	52.829 mL
5 mM	1.0566 mL	5.2829 mL	10.5658 mL
10 mM	0.5283 mL	2.6414 mL	5.2829 mL
50 mM	0.1057 mL	0.5283 mL	1.0566 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481