

CY5.5-Alginate

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	CY5.5-Alginate is a CY5.5-labeled sodium alginate. This compound merges the polysaccharide properties of sodium alginate with the fluorescent characteristics of CY5.5, making it significant in areas such as fluorescence labeling, cell imaging, and drug delivery research (Ex/Em = 680/710 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481