

## CRBN ligand-255

## Chemical Properties

CAS No. :	2925097-84-1
Formula:	C13H12N4O2
Molecular Weight:	256.27
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	CRBN ligand-255 is an E3 ubiquitin ligase cereblon (CRBN) ligand, utilized for the recruitment of cereblon protein. By connecting through a linker, CRBN ligand-255 can bind to a target protein ligand to form a PROTAC.
Targets(IC50)	Ligands for E3 Ligase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9021 mL	19.5107 mL	39.0213 mL
5 mM	0.7804 mL	3.9021 mL	7.8043 mL
10 mM	0.3902 mL	1.9511 mL	3.9021 mL
50 mM	0.078 mL	0.3902 mL	0.7804 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481