

CRBN ligand-435

Chemical Properties

CAS No. :	2925087-53-0
Formula:	C13H11N3O3
Molecular Weight:	257.25
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	CRBN ligand-435 is an E3 ubiquitin ligase cereblon (CRBN) ligand utilized for recruiting the cereblon protein. CRBN ligand-435 can be linked via a linker to a target protein ligand to form a PROTAC.
Targets(IC50)	Ligands for E3 Ligase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8873 mL	19.4363 mL	38.8727 mL
5 mM	0.7775 mL	3.8873 mL	7.7745 mL
10 mM	0.3887 mL	1.9436 mL	3.8873 mL
50 mM	0.0777 mL	0.3887 mL	0.7775 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481