

But-3-yn-1-yl methanesulfonate

Chemical Properties

CAS No. :	72486-09-0
Formula:	C ₅ H ₈ O ₃ S
Molecular Weight:	148.18
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	But-3-yn-1-yl methanesulfonate is a PROTAC linker used in the synthesis of PROTAC molecules.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.7485 mL	33.7427 mL	67.4855 mL
5 mM	1.3497 mL	6.7485 mL	13.4971 mL
10 mM	0.6749 mL	3.3743 mL	6.7485 mL
50 mM	0.135 mL	0.6749 mL	1.3497 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481