

(1S,3S)-3-Benzoic acid 3-carboxycyclopentyl ester

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(1S,3S)-3-Benzoic acid 3-carboxycyclopentyl ester ((1S,3S)-3-(Benzoyloxy)cyclopentanecarboxylic acid) is a synthetic intermediate used in pharmaceutical synthesis.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481