

(+)-Ketodoxapram**Chemical Properties**

CAS No. :	1415394-63-6
Formula:	C ₂₄ H ₂₈ N ₂ O ₃
Molecular Weight:	392.5
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	(+)-Ketodoxapram ((R)-AHR-5955) is an orally active respiratory stimulant. It is applicable in studies related to respiratory system disorders.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5478 mL	12.7389 mL	25.4777 mL
5 mM	0.5096 mL	2.5478 mL	5.0955 mL
10 mM	0.2548 mL	1.2739 mL	2.5478 mL
50 mM	0.051 mL	0.2548 mL	0.5096 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481