

VLVLDTDYK

Chemical Properties

CAS No. : 222406-95-3
Formula: C49H80N10O16
Molecular Weight: 1065.23
Keep away from moisture
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	VLVLDTDYK is a peptide derived from β -lactoglobulin.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9388 mL	4.6938 mL	9.3876 mL
5 mM	0.1878 mL	0.9388 mL	1.8775 mL
10 mM	0.0939 mL	0.4694 mL	0.9388 mL
50 mM	0.0188 mL	0.0939 mL	0.1878 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481