

5-(Benzyloxy)pentan-1-ol

Chemical Properties

CAS No. :	4541-15-5
Formula:	C ₁₂ H ₁₈ O ₂
Molecular Weight:	194.27
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	5-(Benzyloxy)pentan-1-ol is a PROTAC linker utilized in the synthesis of PROTAC molecules.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.1475 mL	25.7374 mL	51.4748 mL
5 mM	1.0295 mL	5.1475 mL	10.295 mL
10 mM	0.5147 mL	2.5737 mL	5.1475 mL
50 mM	0.1029 mL	0.5147 mL	1.0295 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481