

Bid BH3-Gly-r8

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Bid BH3-Gly-r8 is a variant of Bid BH3 linked with eight arginines. It shows no difference in apoptotic effects compared to the unmodified peptide homologue, Bid BH3. Both Bid BH3-Gly-r8 and Bid BH3 exhibit high affinity for Bax and Bcl-2. This compound is applicable in tumor research.
Targets(IC50)	Apoptosis,Bcl-2 Family

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481