

CRBN ligand-471

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | CRBN ligand-471 is an E3 ubiquitin ligase cereblon (CRBN) ligand designed to recruit cereblon proteins. By using a linker, CRBN ligand-471 can be attached to a target protein ligand to form a PROTAC. |
| Targets(IC50) | Ligands for E3 Ligase |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481