

## Polyestradiol phosphate

## Chemical Properties

CAS No. : 28014-46-2

Formula: (C<sub>18</sub>H<sub>24</sub>O<sub>2</sub>.H<sub>3</sub>O<sub>4</sub>P)<sub>x</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Polyestradiol phosphate (Polyestrol phosphate) is a polymer of Estradiol and Phosphoric acid. It is utilized in research related to prostate cancer.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481