

JAK1/2 ligand 1

Chemical Properties

CAS No. :	1026028-96-5
Formula:	C17H13N3O2
Molecular Weight:	291.31
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	JAK1/2 ligand 1 is a ligand for JAK1/2, and can be utilized in the synthesis of PROTACs, such as [JAK1/2 Ligand-Linker Conjugates 1].
Targets(IC50)	JAK,Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4328 mL	17.1638 mL	34.3277 mL
5 mM	0.6866 mL	3.4328 mL	6.8655 mL
10 mM	0.3433 mL	1.7164 mL	3.4328 mL
50 mM	0.0687 mL	0.3433 mL	0.6866 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481