

Desamino lenalidomide-5-F

Chemical Properties

CAS No. :	2468780-76-7
Formula:	C13H11FN2O3
Molecular Weight:	262.236
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Desamino lenalidomide-5-F is an E3 ubiquitin ligase cereblon (CRBN) ligand used to recruit cereblon proteins. It can be linked to a target protein ligand via a linker to form a PROTAC.
Targets(IC50)	Ligands for E3 Ligase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8133 mL	19.0665 mL	38.133 mL
5 mM	0.7627 mL	3.8133 mL	7.6266 mL
10 mM	0.3813 mL	1.9067 mL	3.8133 mL
50 mM	0.0763 mL	0.3813 mL	0.7627 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481