

tert-Butyl (piperidin-4-ylmethyl)carbamate hydrochloride

Chemical Properties

CAS No. :	1049727-98-1
Formula:	C ₁₁ H ₂₃ ClN ₂ O ₂
Molecular Weight:	250.77
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	tert-Butyl (piperidin-4-ylmethyl)carbamate hydrochloride is a PROTAC linker employed in the synthesis of PROTAC molecules.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9877 mL	19.9386 mL	39.8772 mL
5 mM	0.7975 mL	3.9877 mL	7.9754 mL
10 mM	0.3988 mL	1.9939 mL	3.9877 mL
50 mM	0.0798 mL	0.3988 mL	0.7975 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481