

Rhodamine B-Dextran chloride (MW 500000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:
Keep away from direct sunlight
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Rhodamine B-Dextran (MW 500000) is a Rhodamine B-labeled dextran with a molecular weight of 500000. This compound, Rhodamine B-Chitosan (MW 500000), plays a significant role in areas such as fluorescence labeling, cell imaging, vascular imaging, and drug delivery research (Ex/Em = 546/610 nm). |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481