

CRBN ligand-788

Chemical Properties

CAS No. :	2304753-45-3
Formula:	C13H14N4O3
Molecular Weight:	274.28
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	CRBN ligand-788 is a ligand for the E3 ubiquitin ligase cereblon (CRBN), utilized for recruiting the cereblon protein. It can be linked via a linker to a target protein ligand to form a PROTAC.
Targets(IC50)	Ligands for E3 Ligase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6459 mL	18.2295 mL	36.4591 mL
5 mM	0.7292 mL	3.6459 mL	7.2918 mL
10 mM	0.3646 mL	1.823 mL	3.6459 mL
50 mM	0.0729 mL	0.3646 mL	0.7292 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481