

3-(4-Aminophenyl)propionic acid

Chemical Properties

CAS No. :	2393-17-1
Formula:	C9H11NO2
Molecular Weight:	165.189
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	3-(4-Aminophenyl)propionic acid is a PROTAC linker utilized in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.0536 mL	30.2682 mL	60.5364 mL
5 mM	1.2107 mL	6.0536 mL	12.1073 mL
10 mM	0.6054 mL	3.0268 mL	6.0536 mL
50 mM	0.1211 mL	0.6054 mL	1.2107 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481