

Benzyl 5-Bromopentanoate

Chemical Properties

| | |
|-------------------|--|
| CAS No. : | 60343-28-4 |
| Formula: | C ₁₂ H ₁₅ BrO ₂ |
| Molecular Weight: | 271.15 |
| Storage: | Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small> |

Biological Description

| | |
|-------------|--|
| Description | Benzyl 5-Bromopentanoate is a PROTAC linker used in the synthesis of PROTAC molecules. |
|-------------|--|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|-----------|
| 1 mM | 3.688 mL | 18.440 mL | 36.880 mL |
| 5 mM | 0.7376 mL | 3.688 mL | 7.376 mL |
| 10 mM | 0.3688 mL | 1.844 mL | 3.688 mL |
| 50 mM | 0.0738 mL | 0.3688 mL | 0.7376 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481