

## Rhodamine B-Dextran chloride (MW 5000)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
**Store at -20°C**  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Rhodamine B-Dextran (MW 5000) is a Rhodamine B-labeled dextran with a molecular weight of 5000. Rhodamine B-Chitosan (MW 5000) plays a significant role in fluorescence labeling, cell imaging, vascular imaging, and drug delivery research, with excitation/emission wavelengths of 546/610 nm.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481