

N,N'-Dimethyl-1,3-propylenediamine

Chemical Properties

CAS No. :	111-33-1
Formula:	C ₅ H ₁₄ N ₂
Molecular Weight:	102.178
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	N,N'-Dimethyl-1,3-propylenediamine is a PROTAC linker employed in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	9.7867 mL	48.9333 mL	97.8665 mL
5 mM	1.9573 mL	9.7867 mL	19.5733 mL
10 mM	0.9787 mL	4.8933 mL	9.7867 mL
50 mM	0.1957 mL	0.9787 mL	1.9573 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481