

CRBN ligand-203

Chemical Properties

CAS No. :	3023408-98-9
Formula:	C19H22N2O4
Molecular Weight:	342.4
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	CRBN ligand-203 is a ligand for the E3 ubiquitin ligase cereblon (CRBN), used for recruiting the cereblon protein. CRBN ligand-203 can be connected to a target protein ligand via a linker to form a PROTAC.
Targets(IC50)	Ligands for E3 Ligase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9206 mL	14.6028 mL	29.2056 mL
5 mM	0.5841 mL	2.9206 mL	5.8411 mL
10 mM	0.2921 mL	1.4603 mL	2.9206 mL
50 mM	0.0584 mL	0.2921 mL	0.5841 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481