

## Prostaglandin E2-PEG11-biotin

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from direct sunlight

Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Prostaglandin E2-PEG11-biotin is a ligand that can be immobilized on flow cells. Once immobilized, Prostaglandin E2-PEG11-biotin binds BtOCT1. It is applicable for use in surface plasmon resonance (SPR) techniques.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481