

Rhodamine B-Chitosan (MW 3000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Rhodamine B-Chitosan (MW 3000) is a chitosan labeled with Rhodamine B, having a molecular weight of 3000. This compound merges the polysaccharide properties of chitosan with the fluorescent traits of Rhodamine B, playing a significant role in fluorescence labeling, cell imaging, and drug delivery research (Ex/Em = 546/610 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481