

5'-O-DMTr-2'-O-Me-rA

Chemical Properties

CAS No. : 103285-23-0
Formula: C₃₂H₃₃N₅O₆
Molecular Weight: 583.65
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	5'-O-DMTr-2'-O-Me-rA is a monomer used in nucleic acid synthesis.
Targets(IC50)	Nucleoside Antimetabolite/Analog

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7134 mL	8.5668 mL	17.1336 mL
5 mM	0.3427 mL	1.7134 mL	3.4267 mL
10 mM	0.1713 mL	0.8567 mL	1.7134 mL
50 mM	0.0343 mL	0.1713 mL	0.3427 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481